

Oxford Lab Products
BenchMate® Digital Vortex Mixer
 Features and Benefits

Powerful and Broad Speed Range

Adjustable speed range from 300 rpm to 4200 rpm

Easy and Accurate Setting Adjustment

Digital display ideal for protocol documentation

Hands-free Usability

Touch activation and continuous modes for easy mixing

No “Walking”


Built-in counterbalance

Pulse Function

Alternates 5 seconds of mixing with 5 motion-free seconds

Lifetime Warranty*


Durable equipment built to last



Microplate Attachment
(VM-D-PLATE)



24 x 1.5 mL Microtube
Attachment
(VM-D-TUBE)



Round Disk Attachment
(VM-D-ROUND)



Erlenmeyer Flask
Attachment
(VM-D-FLtASK)

Product Specifications for BenchMate® Digital Vortex Mixer	
Catalog Number	VM-D
Speed Range	300 rpm to 4200 rpm
Modes	Touch, Continuous, and Pulse Mode
Timer Function	1 Minute to 999 Minutes and Infinite Mode
Digital Display	Yes
Motor Type	Brushless DC (BLDC) Motor
Orbital Diameter	4 mm
Max. Load Capacity	500 g

*Lifetime warranty only available in the Continental USA