

Oxford Lab Products  
**BenchMate® Mini Vortex Mixer**  
Features and Benefits

**Powerful Vortex Mixing**

2500 rpm to 4500 rpm adjustable speed

**Small Footprint**

Only 100 mm diameter

**Stable with No “Walking” Issues**

Built-in counterbalance and suction cup feet

**Easy to Clean**

Silicone cup

**Robust Design**

Constructed with a die-cast aluminum base ensuring a long-life

**Lifetime Warranty\***

Durable equipment built to last



**Product Specifications for BenchMate® Mini Vortex Mixer**

Catalog Number	VM-M
Speed Range	2500 rpm to 4500 rpm
Modes	Touch Activation
Maximum Tube Diameter	Up to 30 mm
Motor Type	Brushless DC (BLDC) Motor
Orbital Diameter	4.5 mm
Maximum Load Capacity	500 g

\*Lifetime warranty only available in the Continental USA

**Complementary Oxford LP Products**



BenchMate® Centrifuge  
Tubes



BenchMate® Microcentrifuge  
(C8)



Digital Magnetic Stirrer  
(MS1)